


# THE CHATA CLOSED SYSTEM

## THE CHEM+NECT® BAG FOR 100L AND 200L DRUMS


CHATA's CHEM+NECT® products supply the consistency, stability and safety you need for your solutions and provide the increased throughput that helps your lab or manufacturing group reach its full potential.

### THE CHEM+NECT® BAG IS A CLOSED SYSTEM.

The system is constructed of an irradiated multi-layered plastic bag, your sterile filtered solution, and sterilized appurtenant tubing and connectors. All of this placed in a sturdy drum, which may be overwrapped for use in clean rooms or manufacturing suites.



The CHEM+NECT® system prevents exposure of the solution to the environment after manufacture and throughout use. Solutions are removed using a standard peristaltic pump fitted with sterilized tubing that mates easily and effectively with the CHATA CHEM+NECT® bag. As the solution is removed, the bag simply collapses. There is no head space, no exposure to particulates, and no pathway for microorganisms to enter the media and your process. The solution integrity remains unchanged as the drum is being emptied.



CHATA manufactures in ISO certified cleanrooms following ISO and cGMP standards. All solutions and raw materials are fully traceable. The expiry presented on the label is a true expiration date based on actual data, not a warranty.

CHATA'S SOLUTIONS ARE  
CONSISTENT, COMPLIANT, RELIABLE, AND EASY TO USE.